

LogTag Analyser Software and Docking Station Starter Kit



Simple data retrieval from the LogTag recorder. Extensive data analysis and reporting.

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Manufacturer

Description

LogTag Analyser software provides an easy to use powerful platform for configuring the LogTag temperature logger before deployment and for data analysis when data is retrieved. The kit includes the software and the USB docking station.

Logtag Analyser Features:

MultChart Graphing Display.

Hibernate function.

Automatic download of data when a new LogTag is inserted into an interface.

Supports multiple interfaces thereby allow multiple downloads and batch configuring of LogTags.

Display & printing of data in Chart, list or summary form.

Saving downloaded LogTag data to file, retrieving previously filed LogTag data for review and email forwarding of data.

Docking station available with USB(LTI-USB) or RS232(LTI-RS232) serial connection.

Support for Digital Signatures specifically in compliance to FDA 21 CFR Part 11.