

Bench Top Incubator



A small multi-purpose incubator for the modern laboratory, specially for incubation of media and for Compact Dry.

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Manufacturer

Description

Ideal bench top incubator for any laboratory.

Features:

- Holds at least 60 Compact Dry™ bacterial test plates
- Variable Control
- Visual control guaranteed (clear Plexiglass - door)
- Easy to clean
- Includes an adjustable shelf and thermometer
- Can fit an extra shelf

Specifications:

- External Dimensions: 310 x 155 x 168mm (W x H x D)
- Internal Dimensions: 220 x 120 x 150mm (W x H x D)
- Weight (net): 1.1kg
- Consumption: 26W / 220V / 50Hz
- Temperature range: 25 - 45°C (the temperature is factory set at 35°)
- Temperature accuracy: $\pm 1^\circ\text{C}$
- Insulation: Double; SEV tested
- Construction: ABS plastic; white

Ordering Details:

- Product Code: ALINCS70700