

EnviroMax Plus® Environmental Swab



Ideal environmental swab for the microbiological analysis of surfaces and equipment.

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Description

EnviroMax Plus® is a new large oversized sponge swab system, designed to maximise microbial recovery from floors, drains, walls and large machinery surfaces. The swab is pre-moistened with neutralising buffer and peptone water, neutralizing any residual sanitiser that may be present on the test surface.

Pathogen testing is easy using the EnviroMax Plus® swab, as the primary enrichment broth can be simply added to the tube after sampling has been completed.

Advantages

- Quick-Turn-Open Leak Resistant Cap
- Oversized sponge tip for maximum collection
- Pointed sponge tip enables swabbing of corners
- Pre-moistened with a sanitiser neutralising buffer
- Ideal for in-tube primary pathogen enrichment
- Rigid foam paddle tip
- Tip Width/Length – 15.5mm/56mm
- Total Swab Length – 150mm



[Virtual Demonstration - Floor Drain](#)

[Virtual Demonstration - Meat Carcass Sample](#)

[Virtual Demonstration - Stainless Steel Surface](#)

Shelf Life

12 Months from date of manufacture

Storage Conditions

Room Temperature

Ordering Details

Pack of 25

Product Code: PP25-88050PF-50

Pack of 400

Product Code: PP25-88050PFB