

Micro-Touch® Nitra-Tex® Sterile Gloves



Carry out aseptic technique with confidence using the Micro-Touch® Nitra-Tex Sterile Gloves.

Rating: Not Rated Yet
Price

[Ask a question about this product](#)

Description

Micro-Touch® Nitra-Tex® Sterile Gloves provide the user an excellent sensitivity, durability and grip under both wet and dry conditions. This product is made from 100% nitrile and does not contain any natural rubber latex thus addressing any allergy concerns.

Benefits:



- Latex free
- Powder free
- Ambidextrous design
- Textured grip
- Packed in pairs

Quality Standards:

- Meets or exceeds ASTM examination glove standards. Manufactured within the quality guidelines of ISO 13485:2003, ISO 13485:2012, ISO 9001:2000 and FDA-QSR.

Ordering Details:

- Small size: 6034151
- Medium size: 6034152
- Large size: 6034153
- Extra Large size: 6034154