

Warrior® Washdown Scales



Warrior® Washdown Scales is a functional balance for the requirements of the food, manufacturing, pharmaceutical and veterinary industries.

Rating: Not Rated Yet
Price

[Ask a question about this product](#)

Description

Warrior® Washdown Scales is IP-66 rated meeting the demands of washdown and cleanroom environments. The solid stainless steel construction of the base and the indicator provide rugged durability.

Applications:

- Weighing
- Parts counting
- Percentage weighing
- Checkweighing
- Dynamic / animal weighing

Features:

- IP66 rated for wash-down and dust protection
- Vivid, backlit LCD easily visible in any lighting conditions
- Capacity tracker built into display for easily monitoring pending overloads
- Rugged stainless steel (grade 304) construction for ultimate durability
- Rechargeable battery included for operation virtually anywhere
- Programmable auto power-off to save energy
- Selectable digital filtering for animal / dynamic weighing delivers stable results for moving subjects
- RS-232 interface for communication with printers and computers

Capacity (kg)	W500-8	W500-16	W500-32	W500-32	W500-32	W500-75
Capacity	8kg	16kg	32kg	32kg	32kg	75kg
Readability	0.5g	1g	2g	2g	1g	5g
Tare Range	Full range			Full range		
Pan Size	250 x 250mm			300 x 400mm		
Weighing Units	g, kg, lb, oz, N					
Calibration	External calibration					
Display	Backlit LCD with 40mm-high digits					
Power Supply	12VDC 800mA adapter / Internal rechargeable battery					
Operating Temperature	-10°C to 40°C					
Overall Dimensions	250 x 345 x 532mm			300 x 500 x 587mm		

Capacity (kg)	W500-75	W500-150	W500-150	W500-150
Capacity	75kg	75kg	150kg	150kg
Readability	5g	1g	0.01kg	0.01kg
Tare Range	Full range			
Pan Size	400 x 500mm		450 x 600mm	
Weighing Units	g, kg, lb, lb-oz, N			
Calibration	External calibration			
Display	Backlit LCD with 40mm-high digits			
Power Supply	12VDC 800mA adapter / Internal rechargeable battery			
Operating Temperature	-10°C to 40°C			
Operating Temperature	400 x 590 x 972mm		450 x 700 x 972mm	

