

## pH Test Strips



pH Test Strips for the easy and reliable measurements.

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

### Description

pH-Fix test strips are available in a variety of ranges to suit numerous applications across all industries. pH-Fix test strips are supplied for immediate use along with a colour chart for easy result interpretation.

The test strips have four different pads on each strips to cover the pH range. The colour chart matches both the colour and position of the pads on the strips, allowing a precise pH determination.

pH-Fix range	pH gradation	REF
pH-Fix 0 - 14	1 · 2 · 3 · 4 · 5 · 6 · 7 · 8 · 9 · 10 · 11 · 12 · 13 · 14	92110
pH-Fix 0.0 - 6.0	0 · 0.5 · 1.0 · 1.5 · 2.0 · 2.5 · 3.0 · 3.5 · 4.0 · 4.5 · 5.0 · 5.5 · 6.0	92115
pH-Fix 2.0 - 9.0	2.0 · 2.5 · 3.0 · 3.5 · 4.0 · 4.5 · 5.0 · 5.5 · 6.0 · 6.5 · 7.0 · 7.5 · 8.0 · 8.5 · 9.0	92118
pH-Fix 4.5 - 10.0	4.5 · 5.0 · 5.5 · 6.0 · 6.5 · 7.0 · 7.5 · 8.0 · 8.5 · 9.0 · 9.5 · 10.0	92120 CE*
pH-Fix 6.0 - 10.0	6.0 · 6.4 · 6.7 · 7.0 · 7.3 · 7.6 · 7.9 · 8.2 · 8.4 · 8.6 · 8.8 · 9.1 · 9.5 · 10.0	92122
pH-Fix 7.0 - 14.0	7.0 · 7.5 · 8.0 · 8.5 · 9.0 · 9.5 · 10.0 · 10.5 · 11.0 · 11.5 · 12.0 · 12.5 · 13.0 · 13.5 · 14.0	92125
pH-Fix 0.3 - 2.3	0.3 · 0.7 · 1.0 · 1.3 · 1.6 · 1.9 · 2.3	92180
pH-Fix 1.7 - 3.8	1.7 · 2.0 · 2.3 · 2.6 · 2.9 · 3.2 · 3.5 · 3.8	92190
pH-Fix 3.1 - 8.3	3.1 · 3.5 · 3.9 · 4.3 · 4.7 · 5.1 · 5.5 · 5.9 · 6.3 · 6.7 · 7.1 · 7.5 · 7.9 · 8.3	92135 CE*
pH-Fix 3.6 - 6.1	3.6 · 4.1 · 4.4 · 4.7 · 5.0 · 5.3 · 5.6 · 6.1	92130 CE*
pH-Fix 5.1 - 7.2	5.1 · 5.4 · 5.7 · 6.0 · 6.3 · 6.6 · 6.9 · 7.2	92140
pH-Fix 6.0 - 7.7	6.0 · 6.4 · 6.7 · 7.0 · 7.3 · 7.7	92150
pH-Fix 7.5 - 9.5	7.5 · 7.9 · 8.2 · 8.4 · 8.6 · 8.8 · 9.1 · 9.5	92160
pH-Fix 7.9 - 9.8	7.9 · 8.3 · 8.6 · 8.9 · 9.1 · 9.4 · 9.8	92170

\*CE-marked according to the IVQ-directive 98/79/EC



