

24oz (710mL) Whirl-Pak® Stand-Up Bags



24oz Whirl-Pak® Stand-Up Bags provide a convenient means to aseptically take samples.

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Description

Sampling is easy when using the 24oz Whirl-Pak® Stand-Up bag. Once sample is added to the bag, the special gusseted bottom flattens out providing the bag with the ability to stand without any support allowing the user to use both hands.

A white strip is located on the front of the bag, providing an area for sample identification. Ordinary ballpoint pens, pencils and felt tip makers can be used to document the sample identification on the white strip.

Features:

- 'Puncture Proof Tabs' - Tape is extended past the aluminium wire ends eliminating the chances of sharp points puncturing the bag.
- 'Sterile' - Sterilised as a finished product via Ethylene Oxide Gas.
- 'Easy Identification of Sample' - Felt tip and ball point pens, along with pencils can be used to write on bags.
- 'Leak Proof Closure' - By twirling or lightly folding the tab over three times, the Whirl-Pak® bag becomes a leak proof container.

Specifications:

- 24oz (710mL) Capacity
- 15cm x 23cm
- 0.076mm thick
- Box of 500

Capacity is measured in fluid ounces when the bag is closed and tab is folded over three times. Volumes and dimensions are approximate, bags should not be used at temperatures above 82°C. Bags can be frozen to any temperature, even liquid nitrogen - careful handling required after freezing. All bags are sterilised after manufacturing.

Ordering Details:

- Product Code: B01401WA Box of 500