

92oz. Whirl-Pak® Bags



92oz Whirl-Pak® to accommodate for large sample volumes.

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Description

Whirl-Pak® Bags enable easy sterile sampling of liquids, semi-solids and solids for laboratory, product content or legal compliance. This 92oz Whirl-Pak® bag accommodates large size samples and eliminate reducing the need to reduce the sample to smaller sizes.

The closing tab has one round and one flat wire to ensure the bag is closed securely. Fold tab over tightly 4-5 times when closing.

Features:

- 'Puncture Proof Tabs' - tape is extended past the aluminium wire ends eliminating the chances of sharp points puncturing the bag.
- Sterile - Sterilised as a finished product via Ethylene Oxide Gas.
- Easy Identification of Sample - felt tip pen can be used to write on bags.
- Leak proof Closure

Specifications:

- 2721mL (92 ounce) Capacity
- 25.4cm x 38cm
- 0.102cm thick
- Box of 250

Capacity is measured in fluid ounces when the bag is closed and tab is folded over. Volumes and dimensions are approximate, bags should not be used at temperatures above 82°C. Bags can be frozen to any temperature, even liquid nitrogen - careful handling required after freezing. All bags are sterilised after manufacturing.

Ordering Details:

- Product Code: B01445WA Box of 250 bags