

## BagRack® Bag Holder



BagRack® Bag Holder provides an easy and clean means to store / incubate blender bags.

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

### Description



BagRack® Bag Holder is manufactured from stainless steel making it easy to clean and autoclave. The BagRack® is stackable and a drip tray is an optional extra.

The BagRack® is available in four models to suit different bag sizes, refer to ordering details below.

### Features:

- Keeps blender bags upright eliminating spillages.
- Stackable
- Autoclavable

### Ordering Details:

## Accessories: BagRack® Bag Holder

---

- |                         |                       |  |
|-------------------------|-----------------------|--|
| ▪ Product Code: 221 010 | Holds 80 & 100mL Bags | Dimensions (wxdxh): 12.5cm x 27.5cmx14cm   |
| ▪ Product Code: 221 040 | Holds 400mL Bags      | Dimensions (wxdxh): 20.5cm x 35.5cm x 20cm |
| ▪ Product Code: 221 350 | Holds 4 x 3500mL Bags | Dimensions (wxdxh): 43cm x 35cm x 37.5cm   |
| ▪ Product Code: 221 352 | Holds 2 x 3500mL Bags | Dimensions (wxdxh): 35.5cm x 20.5cm x 20cm |