

## MiniMix® 100 Blender Range



MiniMix® 100 range of blenders are ideal for small samples. The adjustable paddles provides a quick and accurate homogenizing of the sample.

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

### Description

MiniMix® 100 range of blenders are specifically designed for small samples up to 10g.

Two models are currently available in the range which have very similar features. The MiniMix® W CC® has a window door and is supplied with a security drip tray, whilst the MiniMix® P CC® has a solid plain door and the security drip tray is an optional purchase. The MiniMix® P CC® is pictured below:



### Ideal for:

- Medical and Veterinary Applications (biopsy, organ tissue sampling)
- Molecular extractions (proteins, lipids, DNA)
- Toxin research

### Advantages:

- Silent and resistant
- Easy Control

- Adjustable speed, time and blending strength
- Optimal bacterial extraction in less than 60 seconds
- Easy insertion of bag with the side-by-side paddles stop
- Easy Cleaning:
  - All stainless steel
  - 270° retractable door
  - Click and Clean® removable paddles

### Specifications:

- All stainless steel
- Plain Door (MiniMix® P CC®), or Window Door (MiniMix® W CC®)
- Multifunction digital screen
- Variable speed (3, 6, 9, 12 strokes / sec)
- Variable time (1 sec. - 1h or ?)
- Variable power, adjustable paddles (patented)
- Q-Tight® closing
- Click & Clean® removable paddles (patented)
- Fully opening door (270°)
- Side-by-side paddle stop
- Useful volume: 5-80mL
- Security drip tray is delivered with the MiniMix® W CC®. (Optional extra for the MiniMix® P CC®)

### Guarantee:

- 3 year guarantee (After registration of the Guarantee form)
- Window door life - time guarantee
- Shock absorbers life-time guarantee