

Reference Material for Gluten Tests



Innovative CRM's to support methods used by food laboratories and food manufacturing companies. Can be used for both Quantitative and Qualitative methods.

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Description

Certified Reference Materials (CRM's) for gluten analysis are now available providing laboratories a user friendly means to verify and validate test methods. Manufactured using intact wheat gluten and certified by multiple ELISA methods, gluten reference materials can be used for both qualitative and quantitative applications.

Formulated at 4 levels: 0, 10, 20 and 40mg/kg as gluten with actual method performance values listed on the Certificate of Analysis.

Ordering Details:

- Blank (0ppm) Product Code: FC-848
- 10ppm Product Code: FC-845
- 20ppm Product Code: FC-846
- 40ppm Product Code: FC-847



[Accredited: ISO Guide 34:2009, ISO/IEC 17025:2005, ISO/IEC 17043:2010, ISO 9001:2008](#)

For further information (including a sample of a Certificate of Analysis) or pricing please contact our office or email sales@arrowscientific.com.au